

Stop chasing down equipment before **Fingertip Revision Amputations**

Eliminate the pre-procedure equipment hunt with SwiftKit™

A purpose built, fingertip revision amputation procedure kit. No OR trays, no SPD reprocessing, nothing to chase down.



Right Instruments, Every Time

No missized rongeurs. No freer with the wrong tip. No digging through multiple trays to find what you need. Just high-quality single-use instruments, every time.

Standardized Setup

30+ minutes saved per case. With **SwiftKit™**, everything is at the bedside in minutes. No OR dependency, no SPD calls, no relying on who's working or what they remember.



**Watch SwiftKit™
In Action**

**ARBUTUS
MEDICAL**



One Pack. Zero Hunting.

ARBUTUS
MEDICAL

What's in SwiftKit™?

QuikRev™ Bone Rongeur

Procedure-specific size and shape. Designed for fingertip revision — not adapted from a general hand tray.




QuikLift™ Freer Elevator

Curved on one side and straight on the other. No more hunting for the right type.



Complete Prep & Dressing Kit

Everything you need for prep and dressing in a single pack. Just grab your lidocaine and gloves.



“While fingertip revision amputations in the ED are an easy procedure, the logistics are a nightmare. I’ve spent more time gathering supplies, such as rongeurs and freers from the OR, than actually performing the procedure. This is frustrating for patients and providers alike, delays patient care, and is a common recurring challenge.”

Greg Pereira MD
Duke University Hospital